



Analytical Instruments

BVSA 

TABLE OF CONTENTS

1

DIFFRACTOMETERS

DRON-8/8T	4
DRON-7M	5
Applications for DRON-8/8T and DRON-7M	6
Software and options for DRON-8/8T and DRON-7M	8
Colibri bench top diffractometer	14

2

SPECTROMETERS AND ANALYZERS

BRA-135F Energy dispersive general purpose spectrometer	17
ASE-2 XRF Energy dispersive sulfur analyzer	22
ASW-2 XRF Wavelength dispersive sulfur analyzer	23
ASE-3 XRF Energy dispersive sulfur analyzer	24
ASW-3 XRF Wavelength dispersive sulfur analyzer	25
AR-35 Automated on-line XRF Wavelength Dispersive slurry analyzer	26
ARF-7 Specialized X-ray Wavelength dispersive analyzer	28

3

X-RAY OPTICAL ELEMENTS AND CRYSTALS-MONOCROMATORS

Optical elements and crystals	30
-------------------------------------	----




DIFFRACTOMETERS

DRON - 8/8T X-RAY MULTIFUNCTIONAL DIFFRACTOMETER

DRON-8/8T X-ray general purpose diffractometer with vertical q-q goniometer and horizontal sample plane enables to perform X-ray diffraction analysis of phase composition, structural state and orientation of widerange of crystalline objects with different shape and size.

- High-precision wide-angle vertical goniometer with changeable radius
- Automated alignment of sample plane
- Implementation of various X-ray diffraction techniques
- Flexible design and wide range of options
- Variety of X-ray optical schemes
- Remote control of all devices and systems




Technical Parameter		DRON-8/8T
Goniometer Type		Vertical θ – θ
X-ray Optical Scheme		Bragg-Brentano / Debye-Scherrer / parallel-beam
Radius R, mm		180 – 250
Angular range, deg	2θ	from -10 to 165
	θ^f	from -5 to 165
	θ^d	from -5 to 95
Scanning modes		discrete / continuous
Scanning methods		θ –0, θ^f , θ^d , Ω , 2θ – Ω , ψ , $\sin^2\psi$
Smallest addressable increment, deg		0.0005 / 0.0001
Scanning rate, deg/min		0.1 – 50
Reproducibility, deg		± 0.001 / ± 0.0001
Maximum angular speed, deg/min		600 / 2000
Radiation doze rate, mSv/h		Not more than 1



DRON - 7M X-RAY MULTIFUNCTIONAL DIFFRACTOMETER

DRON-7M X-ray diffractometer is capable to solve a wide variety of tasks for powder diffraction analysis. Independent control of 2θ and θ movements allows research of single crystals.


- Horizontal two-circle 2θ - θ goniometer
- High reliability and user-friendly operation Flexible design and wide range of options
- High automation for setting and measurements



Technical Parameter		DRON-7M
Goniometer Type		Horizontal 2θ - θ
X-ray Optical Scheme		Bragg-Brentano / Debye-Scherrer / parallel-beam
Radius R, mm		200
Angular range, deg	2θ	from -100 to 165
	θ	from -180 to 180
Scanning modes		discrete / continuous
Scanning methods		$0-2\theta$, 2θ , θ , $2\theta-\Omega$
Smallest addressable increment, deg		0.0001
Scanning rate, deg/min		0.1 – 50
Reproducibility, deg		± 0.0025
Maximum angular speed, deg/min		720
Radiation doze rate, mSv/h		-





APPLICATIONS OF DRON-7M AND DRON-8/8T X-RAY DIFFRACTOMETER

Application Fields	Problems	Samples
Mining industry Mineralogy Construction Machinery Energetics Oil and gas industry Chemistry Electronics Criminalistics Forensics Pharmaceuticals Crystallography Nanotechnology Examination of cultural valuables Ecology	- Qualitative and quantitative phase analysis of polycrystalline materials and objects including coatings and thin films. - Determination of crystallinity, crystallite sizes and microstrains of lattice. - Determination of lattice type and dimensions, crystal structure refinement. - Tracing of phase transitions, chemical reactions and thermal deformations of lattice in variable environment (temperature, pressure, humidity, gaseous medium or vacuum).	
Metallurgy Machinery Electronics Technical crystals	- Analysis of preferred orientation of particles in metals and in other polycrystal- line materials. - Determination of linear, planar and volumetric stresses in welded seams, parts and frameworks. - Determination of orientation of single crystals and different articles made of them.	
Micro- and nano-electronics	- Determination of composition, thickness, mismatch and defects of layers in thin films, epitaxial and nano heterostructures. - Quality control of materials for micro- and nanoelectronics.	
- Catalysis - Colloid chemistry - Electronics - Molecular biology - Automotive- and aircraft industry (plastics and polymers) - Protection of main pipelines and cable industry - Packaging industry (nano composites and films)	Determination of shape, size, phase composition, internal structure, orientation and distribution of nanoparticles in surface-active material, emulsions (including in biological mediums), fibers, catalysts, polymers, nanocomposites, liquid crystals and other disperse systems.	



X-RAY OPTICAL SCHEMES

TYPICAL DIFFRACTION PATTERNS


Debye-Sherrer, grazing incidence and parallel-beam geometries.





Determination of crystal orientation



Single crystals in high resolution geometry



By small-angle X-ray scattering and reflectometry





SOFTWARE FOR DRON-7M AND DRON-8/8T X-RAY DIFFRACTOMETERS

Data processing - DrWin

Processing of diffraction pattern or selection

- Background approximation (by polynomial or user curve)
- Separation of Ka-doublets
- Peak search and determination of their angular positions
- Approximation of reflection profiles by pseudo-Voigt function (for the entire array or independently for each peak)
- Calculation of peak heights and their integral intensities
- Calculation of FWHM of reflections





Quantitative phase analysis - Quan

- Overall analysis of multicomponent mixture
- Analysis of n-component system
- Analysis of sample with known mass absorption coefficient
- Method of internal standard
- Method of Reference Intensity Ratios (RIR's)
- Method of additives
- Method of reduction

Calculation of average size of coherent domains and of microstrains - Size&Strain

- Calculation of size of coherent domains and microstrains by the method of second central moments
- Calculation of instrumental line broadening
- Application of absorption correction to the samples with another composition.



Calculation of theoretical diffraction pattern - TheorPattern

- Simulation of diffraction patterns of multicomponent mixtures from structural data
- Account for of instrumental factor
- Account for texture and crystalline size for each component
- Comparison of simulated and measured diffraction patterns
- Integrated package of geometrical crystallography

Auto indexing of Powder Diffraction Pattern - Ind



- Determination of Bravais lattice type
- Choice of unit cell
- Computation of Miller indices for selected lines
- Bar graph of input diffraction pattern



SOFTWARE FOR DRON-7M AND DRON-8/8T X-RAY DIFFRACTOMETERS

Full profile analysis by Rietveld method - Rietveld

- Refinement of crystal structures from X-ray powder diffraction data of single crystalline phases and mixtures
- Calculation of polynomial and physical background
- Independent refinement of U, V, W, X, Y profile for different phases and for different groups of reflections
- Refinement of unit cell parameters, atomic and thermal parameters, occupations of atomic positions for each phase
- Choice of refinement strategy
- Control of Refinement conditions
- Calculation of five R-factors





Residual stress analysis - MacroStress

- Calculation of peak angular position from center of gravity or from peak topapex
- Application of correction matrix
- Calculation of linear, planar and volumetric stresses
- Calculation of stress deviations

High temperature-X-ray diffraction - Thermo

- 3D-imaging of measured data in "diffraction angle - intensity - temperature" co-ordinates
- Calibration of the measured data set by internal or external standard
- Refinement of unit cell parameters of the calibrated data set
- Determination of phase transition points
- Determination of thermal expansion coefficients (TEC) in different directions and thermal deformation tensors
- Building of TEC figures




Qualitative phase analysis and access to the Powder Diffraction File database - Retrieve and Search-Match

- Use of PDF-2/PDF-4 database of International Center of Diffraction Data (ICDD) for qualitative analysis
- Automatic or manual search algorithm
- Creation of user subbases for search facilitation
- Addition of user standards into subbases
- Qualitative phase analysis by different criteria, bases (subbases)
- Analysis of lines matched by angular position and intensity
- Quantitative phase analysis by Reference Intensity Ratios (RIR's) method
- Access to the data base including search by selected criteria

OPTIONS FOR DRON-7M AND DRON-8/8T DIFFRACTOMETERS

Multidrive attachments and sample stages


Single-axis φ -attachment



For DRON-8/8T

Analysis of textures and residual stresses in polycrystalline materials, determination of single crystal orientation, study of phase composition and structural characteristics of powder and bulk objects.


Two-axis $\varphi\chi$ -attachment



For DRON-7 M, DRON-8/8T

Analysis of lattice dimensions and quality of single crystals in different crystallographic directions.

Four-axis xyz φ holder for large samples



For DRON-8/8T

Mapping of phase composition and structural characteristics across sample surface, analysis of layer composition and structure in thin films.

Maximal sample weight		
1kg	0.3kg	50kg
Maximal sample dimensions (θ x h)		
100x20mm	28x10, 15x100mm	300x250mm
Automatic alignment of sample plane, accuracy		
5 μm	-	5 μm
Smallest selectable step of φ -rotation		
0.001 deg.	0.1 deg.	0.001 deg.
Smallest selectable step of χ -inclination		
-	0.001 deg.	-
Range of χ -inclination		
-	from -3 to 70 deg.	-
Range of xy-movement		
-	-	±100mm
Smallest selectable step of xy-movement		
-	-	0.1mm
Scanning modes		
Ω , Ω - φ , 2θ - Ω , ψ , $\sin^2\psi$	Ω , Ω - φ , χ - φ , 2θ - Ω , ψ , $\sin^2\psi$	Ω , Ω - φ , 2θ - Ω , ψ , $\sin^2\psi$

OPTIONS FOR DRON-7M AND DRON-8/8T DIFFRACTOMETERS

Multidrive attachments and sample stages

Multi purpose five-axis $\chi\varphi\chi$ -attachment



For DRON-8/8T

Analysis of textures and residual stresses in polycrystalline materials, determination of single crystal orientation, study of phase composition and structural characteristics of powder and bulk product.

Mapping of phase composition and structural characteristics across sample surface, analysis layer composition and structure films.

Analysis of lattice dimensions and quality of single crystals in different crystallographic directions, mapping of reciprocal space.

Autosampler



For DRON-7 M, DRON-8/8T

Continuous measurement of powder and bulk samples in automatic mode.
Speed of sample rotation: 0.5 or 1 rps.

Stage for cylindrical samples



For DRON-7M, DRON-8/8T

Measurements of cylindrical samples (capillaries) of 0.1-1.0 mm in diameter in Debye-Scherrer geometry (transition mode).

Maximal sample weight

1kg

-

-

Maximal sample dimensions ($\theta \times h$)

100x10mm

28x25mm

θ 0,1-1.0mm, length up to 100mm

Automatic alignment of sample plane, accuracy

5 μ m

5 μ m

-

Smallest selectable step of φ -rotation

0.001 deg.

-

-

Smallest selectable step of χ -inclination

0.001 deg.

-

-

Range of χ -inclination

from -5 to +95 deg.

-

-

Range of xy-movement

± 20 mm

-

-

Smallest selectable step of xy-movement

0.1 mm

-

-

Scanning modes

Ω , Ω - φ , χ - φ , 2θ - Ω , ψ , $\sin^2\psi$

θ - θ (DRON-8/8T); 2θ - θ (DRON-7M)

2θ

OPTIONS FOR DRON-7M AND DRON-8/8T DIFFRACTOMETERS

X-ray optical elements




One-dimensional parabolic mirror for DRON-7, DRON-8/8T. Converts a divergent primary beam to a parallel one, makes it monochromatic and enhances intensity. Application:

- Measurements of samples with uneven surface.
- Small-angle X-ray scattering (SAXS).
- Grazing-incidence X-Ray diffraction (GIXRD).
- X-Ray reflectivity (XRR).



4-bounce channel-cut Ge 220 x 4 monochromator of Bartels type for DRON-8/8T

- Converts to high-resolution geometry.
- Singles out pure monochromatic Ka1 line with the angular resolution of 12 arc. sec.




Polycapillary lenses for DRON-8/8T Focusing lens provides:

- Intensity gain of primary beam 50-100 times.
- beam spot on sample surface is 50-100 μm .
- Microanalysis in different points of sample surface.


Collimating semi-lens forms quasi-parallel beam of $\varnothing 10 \text{ mm}$ to perform:

- Measurement of uneven surfaces in parallel-beam geometry.
- Analysis of texture and residual stress.



Versatile motorized monochromators for DRON-7M, DRON-8/8T

- Plate or channel-cut crystals.
- Any material (Ge, Si, SiO_2 , LiF, graphite etc.)
- Any crystallographic orientation (111, 100, 110 etc.)
- Any radiation (from Mo to Cr)
- Cuts background and beta-line
- Singles out monochromatic Ka1 line
- Remote adjustment of crystal



OPTIONS FOR DRON-7M AND DRON-8/8T DIFFRACTOMETERS

Fast registration system




Fast registration system with linear stripped PSD Mythen2 R 1D и Mythen2 R 1K (Dectris, Switzerland) for DRON-7M, DRON-8/8T.


- Strip width, μm : 50 ± 3 .
- Number of channels: 1280 (2 R 1K) and 640 (2 R 1D).
- Active area, mm^2 : 8×64 (2R 1K) and 8×32 , 4×32 (2R 1D)
- Measurement time is 100 times less.
- Angular resolution is the same as for scintillation counter.
Increase of signal/noise ratio, especially for the weak reflections.
- Increase of sensitivity limit.
- Suppression of X-ray fluorescence background.
- Automatic calculation of strip aperture, goniometer radius and zero angle during alignment and calibration of PSD

Application:

- Measurements of large number of samples in a limited period of time.
- Analysis of residual stress.
- Study of poorly crystallized and quickly decomposed materials.
- Real-time studies of phase transformations and chemical reactions.
- Identification of minor impurities.
- Measurements of small quantities of material.



Non-ambient chambers



HTK-1200N oven-chamber for DRON-7M, DRON-8/8T Operation temperatures: from +25 to +1200 oC Atmospheres: vacuum (10-4 mbar), air, inert gases.



HTK-16N/2000N strip-heater chambers for DRON-8/8T Tungsten (W) heater (in vacuum): from +25 to +2300 oC Platinum (Pt) heater (in vacuum, on air, or in atmosphere of inert gas): from +25 to +1600 oC

TTK-600 Low-temperature chamber for DRON-8/8T Operation temperatures: from -190 to +600 oC Atmospheres: vacuum (10-2 mbar), air, inert gases

CHC+ cryo & humidity chamber for DRON-8/8T Operation temperatures (in vacuum): from -5 to +400 oC Humidity range: 5 – 95% at temperatures from +10 to +60 oC

Vacuum equipment for DRON-8/8T


Application: tracing of phase transitions and chemical reactions, polymorph screening, analysis of thermal deformations of lattice in variable environment (temperature, pressure, humidity, gaseous medium or vacuum).



COLIBRI BENCHTOP X-RAY DIFFRACTOMETER

COLIBRI portable XRD system is a new advanced solution of Bourevestnik JSC for scientific, educational and industrial applications.

The instrument is a convenient and mobile tool for phase identification and structural analysis of various polycrystalline materials.



- Vertical θ – θ goniometer of unique design with horizontal sample position
- Available configuration with Mythen 1D
- linear position-sensitive detector
- Built-in cooling system
- Pre-aligned and ready to use

Goniometer	Vertical θ – θ
X-ray optical scheme	Bragg-Brentano
Goniometer radius, mm	150
Angular 2θ range, degree	from -5 to +160 (basic configuration) from -3 to +140 (with Mythen 1D PSD)
Slew speed, deg/min	1000
Scanning mode	discrete / continuous
Scanning 2θ speed, deg/min	from 0.01 to 100
Smallest selectable 2θ increment, degree	0.005
Angular accuracy of peak position determination, degree	0.02
High voltage power supply, maximal output power, W	600
Power requirement, V/Hz	220/50 single phase
Power consumption, VA	3500
Weight, kg	100


COLIBRI BENCHTOP X-RAY DIFFRACTOMETER

Software for instrument control, data acquisition and processing of XRD patterns

- Access control to customized functions for different users.
- Push-button execution of high-voltage regime and measurement.
- Build-in touch screen for flexible operation.
- Remote operation and measurement control .
- Easy switch between minimum scan time and maximum resolution to get required data quality.
- Automated data processing after measurement.



Crystallographic software suite for analysis of powder diffraction data (optional)



- Phase identification and quantification of mixtures
- Degree of crystallinity estimation Unit cell determination
- Crystallite size & Lattice strain analysis
- Rietveld refinement of crystal structure

Application fields

- Mineralogy and Mining
- Metallurgy and Machinery
- Cements and Refractories
- Chemistry and Catalysis
- Forensics and Expertise
- Medicine and Pharmacy
- Science and Education
- Environmental control






SPECTROMETERS & ANALYZERS

BRA-135F ENERGY DISPERSIVE GENERAL PURPOSE SPECTROMETER



- Wide range of detected chemical elements – 9F - 92U
- High sensitivity due to optimal X-ray optical path
- High resolution of ultra-modern SDD detector
- Fundamental parameter (FP) Method for steel and alloy quantitative analysis
- Built-in control PC
- Patent for X-ray transparent dividing window



Completely satisfies the requirements of radiation safety.

X-ray fluorescence energy dispersive general purpose spectrometer BRA-135F allows simultaneous determination of chemical elements by characteristic energies in the 1 to 30 keV range (where elements from F till U are fitted) over a wide scope of concentrations from hundreds ppb. BRA-135F analyzes solid, powder and liquid samples, thin layer on the surface or precipitated on filters.

Operating principle


The spectrometer operating principle is based on excitation of fluorescence radiation of atoms in the substance being analyzed by radiation coming from the low-power X-ray tube. The fluorescence radiation from the sample gets into the SDD semiconductor detector where quanta of different energies are converted into electrical pulses, which amplitudes are proportional to the energy of absorbed quanta. Pulse-frequency rate with the certain amplitude is proportional to the chemical element concentration in the sample.

Method advantages

X-ray fluorescence analysis (XRF) occupies a leading position among the other methods for determination of the quantitative elemental composition of substances. XRF advantages are as follows: nondestructive measurements, multielement determination, express method, high accuracy of analysis, wide range of measured concentrations, development level of quantitative analysis theory, possibility for quantitative analysis with absence of standard samples.

High efficiency SDD detector

The silicon drift detector (SDD) with ultra- thin entrance window allows to register x-ray radiation in wide energy range on re- tention of energy resolution and response.



Detector scheme

BRA-135F ENERGY DISPERSIVE GENERAL PURPOSE SPECTROMETER



Changer is able to carry solid, liquid or powder samples



Irregular shape samples

Low detection limit

Owing to optimally selected materials and thickness of primary radiation filters, high transparent X-ray optical scheme, a low detection limit can be achieved for all elements to be analyzed. In the range of light elements from 9F to 17Cl low-energy radiation registration becomes possible using vacuum where the optical path of radiation passes.

Large or irregular shape samples




It is possible to measure large-size or odd-shaped samples:

- Large-size minerals and nuggets;
- Industrial articles for analysis for ROHS requirements;
- Metals and alloys incoming control;
- Analysis of liquids in the special cells or on the special filters.

Compact body and functionality

- The housing spectrometer ensure fully radiation-protection;
- Handles for carrying;
- Built-in operational computer (PC);
- LAN port for remote control and archiving of measurement results
- LIMS integration is available;
- Easy report creation;
- Password protection and separation of access rights

Fields of application



BRA-135F ENERGY DISPERSIVE GENERAL PURPOSE SPECTROMETER

Technical Sheet

Range of detected elements	⁹ F - ⁹² U
Limits of detection without preliminary enrichments, % • for elements from Na to Mg • for elements from Al to Cl • for elements from K to U	$n \cdot 10^{-2}$ 0.002 0.0005
Limit of determination at sample preconcentration (depending on chemical element), %	$1.5 \cdot 10^{-5}$
Limit of determination in mid group element (liquid), g/dm ³	$n \cdot 10^3$
Average time of one sample analysis, s	100
Energy resolution on MnKa line at pulses counting rate up to 10^4 s ⁻¹ , not more than, eV	145
Max. voltage of X-ray tube, kV	50
X-ray power, W	10
X-ray tube cooling	by air
Primary X-ray radiation filters, pcs	5
Number of samples installed into sample changer, up to hanger #1 (Ø34 mm samples)	15
changer #2 (Ø34,36,40,44 mm samples)	11
Maximum sample size, mm	Ø200 x 60
Ethernet connection	Yes
Possibility for remote control	Yes
Overall dimensions (LxHxB), mm	700x410x400
Instrument weight, max, kg	65
Power	220 V, 50 Hz
Power consumption, W	500

BRA-135F ENERGY DISPERSIVE GENERAL PURPOSE SPECTROMETER

Methodology description

Oil analysis

For measuring purposes of trace elements Al, Ba, Ca, Cu, Fe, Mn, V, Ni, Pb, Zn, P into oil and petrochemicals appropriate methodology was developed.

Analytic complex consisting of BRA-135F and measuring methodology is capable of carrying out quantitative element analysis of petrochemicals in order to define metal trace elements and can be used to analyze exhausted motor oils of aircraft, machines, special motor vehicle in order to identify deterioration rate of engines and define applicability of technical service for it. Methodology is purchased additionally.

Besides BRA-135F could be used for testing as per ASTM D4294 and ASTM D6481. These test methods cover the measurement of sulfur, Barium, Calcium, Magnesium, Phosphorus, Zinc & Chlorine in hydrocarbons, such as lubricating used oil, diesel, naphtha, kerosene, residuals, lubricating base oils, hydraulic oils, grease, jet fuels, crude oils, gasoline (all unleaded), and other distillates. Additionally, sulfur in other products, such as M-85 and M-100, may be analyzed using this technique.

Detection limits of BRA-135F according to certified methodology (ppm):

P	Al	Mn	Ba	Pb	V	Cu	Ni	Fe	Zn	Ca
100	100	5	50	5	5	5	5	5	5	50

Cement materials analysis

For measuring purposes of mass fraction of Na, Mg, Al, Si, P, S, Cl, K, Ca, Ti, Cr, Mn, Fe, Zn, Rb, Sr into cements and cement production materials (clinkers, raw mixes) suitable methodology was developed.

Analysis methodology includes algorithms of principal components determination using X-ray energy dispersive fluorescent spectrometers BRA-135F and is based on recommendations from GOST 5382-91, GOST R 55410-2013 (ISO 12677:2011).

There is provided remelting method into platinum crucibles (according to GOST R 55410-2013) for samples preparation in the methodology.

BRA-135F ENERGY DISPERSIVE GENERAL PURPOSE SPECTROMETER



Advanced software

BRA-135F software is the optimal combination of high performance and functionality with an intuitive interface and tooltips. Owing to this, performance of routine measurements requires no special trained staff.

Only few simple steps are enough to make a measurement: select the sample position, enter its code, choose the research method and run the analysis.

Excellent methodical support

Bourevestnik JSC has a methodological support of spectrometers, including the development and validation of methods for determining the elemental composition of various materials: oil products, ores, rock, slag, refining products, cement and raw mixes, soil and sediment, water, air.

Qualitative and quantitative analysis

BRA-135F software unveils wide opportunities for on-line receiving of information on the chemical composition of the material to be analyzed. The user chooses the method of substance analysis: qualitative or quantitative.

Registered spectra can be:

- saved in suitably structured archive;
- got out for repeat analysis based on a new calibration characteristics;
- processed as per user's request: added, deducted, KLM-marked;
- scaled.


For operator convenience in operation, the auxiliary utilities were implemented to minimize errors during the analysis. For example, semi-automatic marker of lines allows correct identification of spectral lines of different elements.

Fundamental parameters method


The software implementing fundamental parameters method allows semi-quantitative and quantitative determination of elements within the range from Mg (12) to Pb (82) in solid samples of steels and alloys with composition of 0.1% to 100%.

If there are no standard composition samples, to calibrate the spectrometer and with large list of materials under analysis, the standardless semi-quantitative analysis is implemented, which represents the following dependences:

- integrated spectrometer sensitivity on Z element;
- relative (to the integral sensitivity) intensity of analytical line corrected for background and drift, on sample composition (this dependence was implemented through absorption parameter P).



Spectrum measuring and its processing conditions



Spectrum measuring and its processing conditions

Element	Symbol	Atomic Number (Z)	Relative Sensitivity	Integration Factor	Background Correction	Drift Correction	Sample Composition Dependence
Aluminum	Al	13	0.15	1.0	0.05	0.02	0.1% to 100%
Silicon	Si	14	0.12	1.0	0.05	0.02	0.1% to 100%
Manganese	Mn	25	0.10	1.0	0.05	0.02	0.1% to 100%
Iron	Fe	26	0.08	1.0	0.05	0.02	0.1% to 100%
Chromium	Cr	24	0.07	1.0	0.05	0.02	0.1% to 100%
Nickel	Ni	28	0.06	1.0	0.05	0.02	0.1% to 100%
Copper	Cu	29	0.05	1.0	0.05	0.02	0.1% to 100%
Zinc	Zn	30	0.04	1.0	0.05	0.02	0.1% to 100%
Gallium	Ga	31	0.03	1.0	0.05	0.02	0.1% to 100%
Germanium	Ge	32	0.02	1.0	0.05	0.02	0.1% to 100%
Arsenic	As	33	0.01	1.0	0.05	0.02	0.1% to 100%
Selenium	Se	34	0.01	1.0	0.05	0.02	0.1% to 100%
Bromine	Br	35	0.01	1.0	0.05	0.02	0.1% to 100%
Krypton	Kr	36	0.01	1.0	0.05	0.02	0.1% to 100%
Barium	Ba	56	0.01	1.0	0.05	0.02	0.1% to 100%
Strontium	Sr	38	0.01	1.0	0.05	0.02	0.1% to 100%
Calcium	Ca	20	0.01	1.0	0.05	0.02	0.1% to 100%
Magnesium	Mg	12	0.01	1.0	0.05	0.02	0.1% to 100%


Fundamental parameters method

ASE-2 XRF ENERGY DISPERSIVE SULFUR ANALYZER



The energy dispersive sulfur (EDX) analyzer for determination of sulfur mass fraction in petrochemicals is according to: EN ISO 20847:2004; ASTM D4294; ISO 13032:2012; EN ISO 8754:2003
Fully radiation-protected.

- Range of determined sulfur concentrations – from 5 mg/kg to 5 %
- Measurement process meets ASTM D4294, ISO 20847
- Helium is not required
- Connection to PC
- LIMS integration




X-ray fluorescence energy dispersive sulfur analyzer ASE-2 is used for the measurement of mass concentration of the sulfur in unleaded gasoline, diesel oil, crude oil, kerosene, petroleum residues, lubricating oil, hydraulic oil, jet engine fuel and other types of cutter oil. X-ray radiation of low-powered X-ray tube converted by primary radiation filter excites atoms fluorescence radiation of the sample being analyzed. Radiation beams (primary X-ray radiation scattered on the sample and secondary fluorescence one) are fed to the gas proportional counter; in this case the fluorescence radiation of sulfur atoms (SK α) is separated from radiation of other energies with the help of selective filters. Intensity of fluorescence radiation of sulfur atoms registered by the gas proportional counter is proportional to sulphur mass fraction in the sample.

Technical Sheet

Sulfur mass fraction determination method	X-ray fluorescence energy dispersive sulfur (EDX) analyzer with selective filters.
Statistic limit of detection, max., ppm	3
Range of determined Sulphur concentrations, ppm	5-50000
Limits of basic relative error, %	± 0.5
Power consumption, W(220 ACV, 50 Hz Mains)	60
Instrument weight, max kg	12

ASE-2 XRF WAVELENGTH DISPERSIVE SULFUR ANALYZER

- Range of determined sulfur concentrations – from 3 mg/kg to 5 %
- Vacuumized measurement chamber, helium purging is not required
- Helium purging option is available
- Touch screen display
- LIMS integration
- Results storage



The wavelength energy dispersive sulfur (WDX) analyzer for determination of sulphur mass fraction in petrochemicals is according to: EN ISO 20884:2004; ASTM D 6334, ASTM D 2622

Fully radiation-protected.

X-ray wavelength dispersive sulfur analyzer ASW-2 is used for the measurement of mass concentration of the sulfur in unleaded gasoline, diesel oil, crude oil, kerosene, petroleum residues, lubricating oil, hydraulic oil, jet engine fuel and other types of cutter oil. Analyzer ASW-2 allows to measure mass concentration of the sulfur in vacuumized measurement chamber mode, as well as in helium purging mode. For this purpose the instrument is equipped with the set of tooling for the connection to helium station.

Technical Sheet

Sulfur mass fraction determination method	X-ray fluorescence wavelength energy dispersive (WDX) analyzer with vacuumized chamber
Statistic limit of detection, max., ppm	1.5
Range of determined Sulphur concentrations, ppm	3-50000
Limits of basic relative error, %	±0.5
Power consumption, W(220 ACV, 50 Hz Mains)	250
Instrument weight, max kg	45


ASE-3 XRF ENERGY DISPERSIVE SULFUR ANALYZER



The energy dispersive sulfur (EDX) analyzer for determination of sulfur mass fraction in petrochemicals is according to: EN ISO 20847:2004; ASTM D4294; ISO 13032:2012; EN ISO 8754:2003

Fully radiation-protected.

- Range of determined sulfur concentrations – from 5 mg/kg to 5 %
- Helium is not required but helium purged optical path to maximize sensitivity is available
- Touch screen display
- Connection to PC
- LIMS integration
- Autosampler



X-ray fluorescence energy dispersive sulfur analyzer ASE-3 is used for the measurement of mass concentration of the sulfur in unleaded gasoline, diesel oil, crude oil, kerosene, petroleum residues, lubricating oil, hydraulic oil, jet engine fuel and other types of cutter oil. X-ray radiation of low-powered X-ray tube converted by primary radiation filter excites atoms fluorescence radiation of the sample being analyzed. Radiation beams (primary X-ray radiation scattered on the sample and secondary fluorescence one) are to get to the gas proportional counter; in this case the fluorescence radiation of sulfur atoms (SK α) is separated from radiation of other energies with the help of selective filters. Intensity of fluorescence radiation of sulfur atoms registered by the gas proportional counter is proportional to sulphur mass fraction in the sample.

Technical Sheet

Sulfur mass fraction determination method	X-ray fluorescence energy dispersive sulfur (EDX) analyzer with selective filters.
Statistic limit of detection, max., ppm	3
Range of determined Sulphur concentrations, ppm	5-50000
Limits of basic relative error, %	± 0.5
Number of samples installed into sample changer, up to	3
Power consumption, W(220 ACV, 50 Hz Mains)	60
Instrument weight, max kg	12


ASE-3 XRF WAVELENGTH DISPERSIVE SULFUR ANALYZER

- Range of determined sulfur concentrations – from 3 mg/kg to 5 %
- Vacuumized measurement chamber, helium purging is not required
- Helium purging option is available
- Touch screen display
- LIMS integration
- Results storage
- Autosampler



The wavelength energy dispersive sulfur (WDX) analyzer for determination of sulphur mass fraction in petrochemicals is according to: EN ISO 20884:2004; ASTM D 6334, ASTM D 2622, ASTM D 4927

Fully radiation-protected.




X-ray wavelength dispersive sulfur analyzer ASW-3 is used for the measurement of mass concentration of the sulfur in unleaded gasoline, diesel oil, crude oil, kerosene, petroleum residues, lubricating oil, hydraulic oil, jet engine fuel and other types of cutter oil. Analyzer ASW-3 allows to measure mass concentration of the sulfur in vacuumized measurement chamber mode, as well as in helium purging mode. For this purpose the instrument is equipped with the set of tooling for the connection to helium station.

Technical Sheet

Sulfur mass fraction determination method	X-ray fluorescence wavelength dispersive sulfur (WDX) analyzer with vacuumized chamber.
Statistic limit of detection, max., ppm	1.5
Range of determined Sulphur concentrations, ppm	3-50000
Limits of basic relative error, %	±0.5
Number of samples installed into sample changer, up to	9
Power consumption, W(220 ACV, 50 Hz Mains)	250
Instrument weight, max kg	45

AR-35 AUTOMATED ON-LINE XRF WAVELENGTH DISPERSIVE SLURRY ANALYZER



The analyzer design provides overall protection of operation personnel against X-ray radiation.

AR-35 analyzer is designed for on-line flow stream X-ray fluorescence analysis of solutions, suspensions and slurries of ore processing. AR-35 analyzer simultaneously measures concentrations of up to 8 chemical elements in technological product that saves time and reduces the cost of element determination.

The analyzer operation principle is based on the excitation of atomic fluorescence radiation of substance sample by X-ray tube radiation. The fluorescence radiation of various chemical elements is dispensed by the analyzing crystal, and then radiation of a particular wavelength is registered by X-ray detector. Intensity of fluorescence radiation registered of a certain wavelength is directly proportional to the mass fraction of a chemical element in the substance tested.

Technical Sheet	
Range of determined chemical element	$^{20}\text{Ca} - ^{92}\text{U}$
Number of simultaneously defined chemical elements	7
Number of flow measuring cells (sequentially analyzed products, flows) per one unit, depending on customer's requirements.	6,12 or 15
Limits of basic relative error, %	$\pm 0.5\%$
Detection limits, ppm - in solutions - in suspensions and slurries	$10 \dots n \cdot 0.1$ $500 \dots 10$
Average time of one flow analysis, s	20-100
Power consumption, W	5
Instrument weight, max, kg	1200

AR-35 AUTOMATED ON-LINE XRF WAVELENGTH DISPERSIVE SLURRY ANALYZER

Distinctive features


- At concentrating mills of mining enterprises of ferrous and non ferrous metallurgy, in chemical industry, etc.
- Compatible with systems of automated sampling, sample delivering, processing and presentation of analysis results, organization of data analysis archive storage.
- Communication with factory automated process control system.
- High efficiency , accuracy of the analysis, low detection limit, analysis reproducibility.
- High reliability

Scope of application

- Automated system of analytical monitoring and SCADA of nonferrous metals ore processing plants (Fe-Cu-Zn-Pb, Fe-Ni-Co-Ni, Cu-Mo, Mo-W) with branched mixed flotation circuits.
- Hydrometallurgical limits of extraction and refining of non-ferrous, rare and scattered elements (Co, Ni In, Ti, Sc, Y, rare earth elements, Nb, Ta, Mo, W, Re, U).

Functional diagram of the analytical control based on AR-35


Automated workstations are based on advanced software and hardware. The software implements two main methods for quantitative analysis: on regression equations and method of standard scattered radiation.



ARF-7 SPECIALIZED XRF WAVELENGTH DISPERSIVE ANALYZER



High spectral resolution at short wavelengths of X-ray radiation (up to 40 eV)



The instrument design provides overall protection of operation personnel against X-ray radiation.

Specialized wavelength energy dispersive analyzer constructed by Cauchois scheme is designed for high- precision quantitative determination of U, Th, Mo, Au, W, Tl, As, Pb, as well as other elements in the ore, rocks and waste deposits development.

Technical Sheet

Range of determined elements	$^{27}\text{Co} - ^{92}\text{U}$ $^{27}\text{Co} - ^{58}\text{Ce}$, K-series radiation, $^{73}\text{Ta} - ^{92}\text{U}$, L-series radiation.
Spectral resolution (half-width of line U $\text{L}\alpha_1$) less than, eV	40
Range of determined concentrations, %	from 10^{-4} to 100
Limits of basic relative error, %	0.5
Detection limit of ^{92}U for 100 s, max, ppm	1.5
Number of simultaneously loaded samples, pcs.	30
Power of consumption, kW	4.6
Instrument weight, kg	400



ARF-7 SPECIALIZED XRF WAVELENGTH DISPERSIVE ANALYZER

The analyzer operation principle is based on the excitation of atomic fluorescence radiation of substance sample by X-ray tube radiation. Fluorescence emission can be decomposed in spectrum by Cauchois method.

Fluorescence radiation focused by the analyzing crystal of the determined element and the line of standard are allocated onto the Rowland circle. Then they are recorded by turns by the X-ray detector. Intensity of fluorescence

radiation registered of a certain wavelength is directly proportional to the mass fraction of a chemical element in the substance tested.


Distinctive scopes of application

1. Geology and mining industry:
 - determination of rock and ore composition;
 - determination of concentrates composition.
2. Environmental safety:
 - determination of Co, Ni, Cu, Zn, Ga, As, Rb, Sr, Ba, Pb in soils, sediments and rocks at ppm contents.
3. Determination of uranium and new advanced nuclear fuel thorium in rocks, ores and technological products in a wide range of concentrations from 10-4% and above.
4. Determination of uranium in rocks, ores and products of their processing by X-ray fluorescence measurements is carried out according to the procedure No.420-PC developed by VIMS.


Distinctive scopes of application

- Ability to define groups of elements without reconfiguring analyzing crystal.
- Extremely high resolution of Cauchois X-ray optical
- scheme with quartz analyzing crystal(1011).
- Intelligent mathematical support.

Program Interface



Calibration diagram for low concentrations of uranium in geological samples.



Calibration diagram for high concentrations of uranium in geological samples.



X-RAY OPTICAL ELEMENTS AND CRYSTALS- MONOCHROMATORS

X-RAY OPTICAL ELEMENTS AND CRYSTALS-MONOCROMATORS

Crystal diffraction dispersing elements (CDEs)

Crystal diffraction dispersing elements (CDEs) are designed for X-ray dispersion in equipment for spectral (analyzing crystals) and structural analysis (monochromator crystal).

“Bourevestnik” JSC manufactures CDEs optimized with respect to diffraction parameters for serial instruments produced by domestic and world’s leading manufacturers and by the orders of research centres.



Asymmetrical Johansson Monochromator, R188, cut angle 10°59', LiF (220), 2d=0.284 nm

Distinctive features

Materials and reflecting planes of CDEs

- lithium fluoride (200), (220), (420)
- germanium (111), (220), (422)
- silicon (111), (220), (422)
- quartz (1011), (1010), (1340) • graphite (0002)
- pentaerythritol PET(002)
- potassium biphthalate KAP(001)
- rubidium biphthalate RbAP(001)
- multilayer coatings of different types
- other materials and reflecting planes at the customer's request

Different orientation of reflecting plane:

- parallel to working surfaces (reflection geometry)
- parallel to working surfaces (transmission geometry)
- an optional angle to the working surface

CDEs of complex configuration:

- single bent (cylindrical, conical)
- double bent (spherical, toroidal)
- with complex contour (triangular, elliptical, etc.)
- cylindrically grinded and bent (Johansson type)
- bent to logarithmic spiral
- channel-cut for multiple reflection and interference



Johann cylindrical analyzer. InSb(111), 2d=0,748 nm



Flat monochromator. Graphite (0002), 2d= 0,668 nm



Toroidal monochromator. R1=150, R2=75, LiF(200), 2d=0,403 nm

X-RAY OPTICAL ELEMENTS AND CRYSTALS-MONOCROMATORS

Single-crystal sample holders

Single-crystal sample holders are designed for diffraction analysis of weak reflections and microquantity of samples on diffractometers of various brands including DRON. Due to interference in the single-crystal, they eliminate coherent scattering which is the main part of the registered background while using other holders.

Distinctive features

- Materials: silicon, quartz
- Shape: round, rectangular
- Working surface: flat polished or with a sample cavity (cell)



*X-ray interferometer.
Si(220), $2d=0,384\text{ nm}$*



*Spherical analyzer.
R500, KAP(001), $2d=2,664\text{ nm}$*



Flat silicon single crystal holders.



Single-crystal silicone holders with a cavity.